

Product Datasheet

Anti-Hu CD89 Purified (orb371925)

Description	Mouse monoclonal antibody to CD89
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Ag8.653 myeloma cells
Target	CD89
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Clone Number	A59
Antibody Type	Primary Antibody
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P24071
Entrez	2204

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Expiration Date

12 months from date of receipt.



Separation of human neutrophil granulocytes (red-filled) from CD89 negative lymphocytes (black-dashed) in flow cytometry analysis (surface staining) of human peripheral whole blood stained using anti-human CD89 (A59) purified antibody (concentration in sample 3 µg/ml) GAM APC.



Flow cytometry surface staining pattern of human peripheral whole blood stained using anti-human CD89 (A59) purified antibody (concentration in sample 3 µg/ml) GAM APC.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)