

Product Datasheet

Creatine Kinase Activity Colorimetric Microplate Assay Kit (orb390751)

Catalog Number	orb390751
Category	Assays and Kits
Description	Detection and Quantification of Creatine Kinase Activity.
Reactivity	All
Tested applications	FA
Application notes	Detection Method: Colorimetric
Assay Type	Detection
Range	30 $\mu\text{mol/L}$ - 600 $\mu\text{mol/L}$
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-15-4
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Representative Standard Curve



Positive Control reaction in 96-well plate assay with decreasing the concentration

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)