

Product Datasheet

Hot Start Polymerase Ab+ Kit (orb533337)

Catalog Number orb533337

Category Assays and Kits

Description Hot Start Polymerase Ab+ provides improved specificity and sensitivity when

room-temperature set up is required. The polymerase activity is blocked at ambient temperature and switched on automatically at the initial denaturation. The ther-mal activation prevents the extension of nonspecifically annealed primers and primer-dimer formation at low temperatures during PCR setup. The polymerase is recommended for routine PCR applications (up to 4 kb fragment

amplifying low-copy-number targets in complex backgrounds or when prolonged

length), high throughput PCR or genotyping. The Crystal Buffer system guarantees robust and reliable amplifi-cation results in almost all PCR applications. The buffer contains a well-balanced ratio of potassium-, ammonium-and magnesium-ions to ensure high specificity and minimal byproduct formation without the need of additional optimization steps. Ruby Buffer

additionally contains gel loading buffer and an inherent red dye allowing the direct loading of the PCR product into the gel. The red dye allows an easy visual control during PCR set-up and in combination with the density reagent the direct loading of the reaction product into the gel. The enzyme replicates DNA at 72 °C. It catalyzes the polymerization of nucleotides into duplex DNA in 5'→3' direction

in the presence of magnesium. It also possesses a $5'\rightarrow 3'$ polymerization-dependent exonuclease replacement activity but lacks a $3'\rightarrow 5'$ exonuclease

(proof-reading) activity.

Form/Appearance liquid

Concentration 5 units/μl

Storage store at -20 °C

Note For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road Cambridge CB5 8LA United Kingdom

Email: <u>info@biorbyt.com</u>, <u>support@biorbyt.com</u> Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive Research Triangle Park Durham NC 27713 United States

Email: <u>info@biorbyt.com</u>, <u>support@biorbyt.com</u> Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558