

## Product Datasheet

### SCRIPT RT-PCR Master (2x) (orb692273)

<b>Catalog Number</b>	orb692273
<b>Category</b>	Tools
<b>Description</b>	<p>SCRIPT RT-PCR Master is designed for performing highly sensitive and specific RT-PCR convenient in single tubes. The enzyme mix is based on a genetically engineered reverse transcriptase with enhanced thermal stability providing increased specificity, high cDNA yield and improved efficiency for highly structured and long cDNA fragments. The 2 x conc. mix contains all reagents required for RT-PCR (except template and primer) in one tube to ensure fast and easy preparation with a minimum of pipetting steps. Premium quality enzymes in combination with ultrapure dNTPs and an advanced buffer system ensure superior amplification results. RT-PCR is used to amplify double-stranded DNA from single-stranded RNA templates. In the RT step the reverse transcriptase synthesizes single-stranded DNA molecules (cDNA) complementary to the RNA template. In the first cycle of the PCR step hot start DNA polymerase synthesizes DNA molecules complementary to the cDNA, thus generating a double-stranded DNA template. During subsequent rounds of cycling the DNA polymerase exponentially amplifies this double-stranded DNA template. In one-step RT-PCR all components of RT and PCR are mixed in one tube so that the complete reaction can be performed without opening the tube. This offers tremendous convenience when applied in routine testing and minimizes the risk of contaminations. SCRIPT RT-PCR Master contains RNase inhibitor that is essential when working with low amounts of starting RNA.</p>
<b>Form/Appearance</b>	liquid
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558