

Product Datasheet

human beta-Endorphin ELISA Kit (orb405929)

Description	human Beta-Endorphin,beta-EP ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Tested Applications	ELISA
Target	β-EP
Note	For research use only.
Sample Types	serum, urine, cerebrospinal fluid (CSF), saliva, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, urine, cerebrospinal fluid (CSF), saliva, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm