

Product Datasheet

Human PSME3 ELISA Kit orb406183

Description	Human Proteasome activator complex subunit 3(PSME3) ELISA kit
Reactivity	Human
Range	18.75 pg/ml-1200 pg/ml
Sensitivity	4.68 pg/ml
Tested Applications	ELISA
Target	PSME3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm