

Product Datasheet

Dog β -EP ELISA Kit orb403512

Description	Canine Beta-Endorphin, β -EP ELISA Kit
Reactivity	Canine
Range	20 pg/ml- 4000 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Target	β -EP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm