
Product Datasheet

Horse β -Endorphin ELISA Kit (orb403848)

Description	Horse Beta-Endorphin, β -EP ELISA Kit
Reactivity	Equine
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.78 pg/mL
Tested Applications	ELISA
Target	β -EP
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm