

Product Datasheet

TUNEL BrightGreen Apoptosis Detection Kit (A112-01)

Description	During apoptosis, intracellular specific endonucleases are activated, chromatin DNA is specifically ...
Storage	Store at -30 ~ -15°C and transport at ≤0°C. BrightGreen Labeling Mix: Store at -30 ~ -15°C and protect from light.
Note	For research use only
Application notes	Bright factor: Brighter fluorescence, stronger anti-fade ability Efficient: Highly active recombinant TdT enzyme ensures fluorescence incorporation efficiency Broad compatible: Suitable for cell slides, paraffin sections and frozen sections
Expiration Date	Please enquire.

