

Product Datasheet

Annexin V-FITC/PI Apoptosis Detection Kit (A211-01)

Description

Annexin V is labeled with green fluorescent FITC, and the labeled Annexin V retains a high affinity ...

Storage

Store at 2 ~ 8°C and protect Annexin V-FITC and PI Staining Solution from light. If not in use for a long time, store Annexin V-FITC at -30 ~ -15°C to avoid repeated freezing and thawing. Adjust the shipping method according to the destination.

Note

For research use only

Application notes

Double label: Annexin V-FITC and PI; distinguish cells at different stages of apoptosis (early / late apoptotic) Additional features: Red fluorescent label for easy co-labeling with other green fluorescents such as GFP High binding efficiency: Annexin V conjugates have higher binding efficiency and extremely accurate characterization of the apoptotic process

Expiration Date

Please enquire.

