

Product Datasheet

CCK-8 Cell Counting Kit (A311-02)

Description	CCK-8 Cell Counting Kit is based on WST-8 [2-(2-methoxy-4-nitrophenyl)-3-(4-nitrophenyl)-5-(2,4-disu...
Storage	Store at 2 ~ 8°C and protect from light. Adjust the shipping method according to the destination.
Note	For research use only
Application notes	Simple operation: Ready-to-use reagent Reliable performance: Low cytotoxicity, high sensitivity, and good repeatability
Expiration Date	Please enquire.

