

Product Datasheet

Vazyme - ClonExpress MultiS One Step Cloning Kit (C113-02)

Description	The ClonExpress technology is simple, fast, and highly efficient DNA seamless cloning technology. It...
Storage	Store at -30 ~ -15°C and transport at ≤0°C. ▲Please avoid repeated freezing and thawing.
Note	For research use only
Application notes	Good performance: No ligase, no self-ligation, positive rate > 95% Broad Compatibility: Efficient cloning of 50 bp - 10 kb DNA fragments; applicable for 2 - 5 fragments Additional features: No need to consider the restriction site of insert fragments Software support: Provide online primer design software "CE Design"
Expiration Date	Please enquire.

