

## Product Datasheet

### Myco-Off Mycoplasma Cleaner (D103-01)

|                          |                                                                                                                                                                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | This product consists of double-stranded DNA fragments of specific molecular weight and is mixed wit...                                                                                                                                                                                      |
| <b>Storage</b>           | Store at -30 ~ -15°C and transport at ≤0°C.                                                                                                                                                                                                                                                  |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                        |
| <b>Application notes</b> | Good stability: No change in product performance at room temperature for 6 months<br>Reliable results: Electrophoretic bands are Sharp and clear<br>Convenient: Provided with Gel Loading Buffer for tracking of sample DNA migration<br>Precise: Exact molecular weight of DNA in each band |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                                                              |

