

Product Datasheet

SuperPico ECL Chemiluminescence Kit (E422-01)

Description Vazyme SuperPico ECL Chemiluminescence Kit is a sensitive, luminol-based enhanced chemiluminescent s...

Storage Store at 2 ~ 8°C and protect from light. Adjust the shipping method according to the destination.

Note For research use only

Application notes High sensitivity: The detection threshold of protein amount reaches picogram level. Strong anti-overexposure ability: The overexposure time of the high microgram protein reaches 15 min. Excellent long-Lasting fluorescence ability: Intact bands can be observed after 4 h.

Expiration Date Please enquire.

