

Product Datasheet

SuperFemto ECL Chemiluminescence Kit (E423-01)

Description Vazyme SuperFemto ECL Chemiluminescence Kit is a sensitive, luminol-based enhanced chemiluminescent ...

Storage Store at 2 ~ 8°C and protect from light. Adjust the shipping method according to the destination.

Note For research use only

Application notes High sensitivity: The detection threshold of protein amount reaches 32 fg. Strong anti-overexposure ability: The overexposure time of the high microgram protein reaches 2 min. Short automatic exposure time: Hundreds of milliseconds. Excellent long-Lasting fluorescence ability: Intact bands can be observed after 4 h.

Expiration Date Please enquire.

