

Product Datasheet

VAHTS RNA Multiplex Oligos Set 2 for Illumina (N324-01)

Description	T4 DNA Ligase catalyzes the formation of phosphodiester bonds between juxtaposed 5' phosphate and 3'...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Quality Control: > 99% purity confirmed by HPLC; nuclease-free components Reliable performance: Classical products for efficient ligation reactions
Expiration Date	12 months from date of receipt.

