

Product Datasheet

VAHTS mRNA Capture Beads 2.0 (N403-01)

Description	VAHTS mRNA Capture Beads 2.0 are 1 µm Oligo(dT)-coupled paramagnetic microspheres intended for the ...
Storage	Store at 2 ~ 8°C. Adjust the shipping method according to the destination.
Note	For research use only
Application notes	High purity: Specific enrichment of mRNA, higher purity and reduced rRNA residue. Wide compatibility: Compatible with multiple species - eukaryotes such as animals, plants, fungi. Compatible with a wider range of input samples—10 ng-12.5 µg.
Expiration Date	12 months from date of receipt.

