

## Product Datasheet

# VAHTS PCR-Free DNA Adapters Set 4 for MG (NM10904)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | VAHTS PCR-Free DNA Adapters for MGI is a specialized kit for MGI high-throughput sequencing platform... |
| <b>Storage</b>           | Store at -30 ~ -15°C and transport at ≤0°C.                                                             |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | Apply to MGI platform for PCR-Free DNA library preparation. Providing 96 different single indexes.      |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                         |

