

Product Datasheet

Library Dilution Buffer (NQ106)

Description VAHTS Library Quantification Kit for Illumina is specially designed for accurate quantification of t...

Storage Store at -30 ~ -15°C and transport at ≤0°C. Master Mix and ROX should be protected from light.

Note For research use only

Application notes Contains all components quantified. New dye method qPCR master mix VAHTS SYBR qPCR Master Mix. High amplification efficiency and good GC content compatibility. Accurate quantification of standard products, CV value between batches 5%. Carefully optimized ROX dye concentration, providing ROX premixed version. Compatible with all mainstream quantitative PCR instruments.

Expiration Date 12 months from date of receipt.

