

Product Datasheet

3G Taq Master Mix for PAGE (Red Dye) (P115-02)

Description	3G Taq Master Mix for PAGE (Red Dye) contains 3G Taq DNA Polymerase, dNTP Mix, visualization red dye...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Eye-catching colour: Red tracking dye, the colour is more striking Convenient operation: Ready-to-use mix; pre-mixed formulation saves time and reduces contamination Robust performance: High amplification efficiencies and inhibitor resistance
Expiration Date	12 months from date of receipt.

