

## Product Datasheet

### 2 × Taq Plus Master Mix (P211-03)

|                          |                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | 2 × Taq Plus Master Mix contains Taq Plus DNA Polymerase, dNTPs, and an optimized buffer system, wh...                                                                                                                                                                                                      |
| <b>Storage</b>           | Store at -30 ~ -15°C and transport at ≤0°C.                                                                                                                                                                                                                                                                 |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                       |
| <b>Application notes</b> | Good performance: ~6X higher than Taq DNA Polymerase; higher amplification yield<br>Convenient operation: Ready-to-use mix; pre-mixed formulation saves time and reduces contamination<br>Additional features: 3'→5' exonuclease activity<br>Good Stability: Stable performance after 50 freeze-thaw cycles |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                                                                             |

