

Product Datasheet

2 × Taq Plus Master Mix II (Dye Plus) (P213-01)

Description	2 × Taq Plus Master Mix II (Dye Plus) contains Taq Plus DNA Polymerase, dNTP and an optimized buffer...
Storage	Store at -30 ~ -15°C for up to 12 months, and transport at ≤0°C.
Note	For research use only
Application notes	Add 3'→5' exonuclease activity; Fidelity is 6 times that of Taq DNA Polymerase; Improve amplification yield and length; Ready-to-use premix, minimizing experimental steps
Expiration Date	12 months from date of receipt.

