

## Product Datasheet

### Vazyme - Plant Direct PCR Kit (PD105-02)

<b>Description</b>	Plant Direct PCR Kit is suitable for direct amplification of plant leaves, seeds, etc., and can be u...
<b>Storage</b>	2 × Plant Direct Master Mix should be stored at -20°C and avoid repeated freezing and thawing; Plant Direct Lysis Buffer can be stored at -20°C or 2°C - 8°C for up to 18 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Special lysate based on proteinase K; The lysate is directly used for PCR amplification without the need to extract the genome; It is widely suitable for the amplification of target fragments within 2 kb, and is suitable for multiple PCR reactions within 4 pairs of primers.
<b>Expiration Date</b>	Please enquire.

