

Product Datasheet

AceQ Universal SYBR qPCR Master Mix (Q511-02)

Description AceQ Universal SYBR qPCR Master Mix is a special premix for qPCR detection using SYBR Green I chimer...

Storage Store at -30 ~ -15°C and protect from light. Transport at ≤0°C.

Note For research use only

Application notes Universal ROX Passive Reference Dye: It is suitable for all qPCR instruments. The concentration of ROX does not need to be adjusted on different instruments. Rigorous hot-start DNA polymerase: Novel chemically modified hot-start Ace Taq DNA polymerase with strong amplification performance and high specificity. Optimized reaction system: Inhibit non-specific amplification and primer dimerization to a large extent. Stable well-to-well repeatability: Meet the requirements of high-class magazines.

Expiration Date 12 months from date of receipt.

