

Product Datasheet

AccurSTART One Step RT-qPCR Probe Kit (Glycerol-free) (QL227-EN02)

Description AccurSTART One Step RT-qPCR Probe Kit (Glycerol-free) is a glycerol-free reagent using RNA as a temp...

Storage Store at -30 ~ -15°C and transport at ≤ 0°C.

Note For research use only

Application notes Sensitive and efficient: Excellent amplification performance
Excellent tolerance to impurities: Compatible with different samples (e.g. feces, tissue, and swab)
Strong lyophilizable stability: Superior stability of accelerated test and freeze-thaw cycles

Expiration Date Please enquire.

