

Product Datasheet

Taq Pro Multiple Probe qPCR Mix (QN213-EN02)

Description Taq Pro Multiple Probe qPCR Mix is a master mix for probe qPCR to detect DNA templates (such as DNA ...

Storage Store at -30 ~ -15°C and transport at ≤0°C.

Note For research use only

Application notes Great low-copy templates amplification sensitivity and 'S-shaped' amplification curve
Compatible with Fast Program Compatible with multiplex detection Strong anti-contamination ability

Expiration Date 12 months from date of receipt.

