

Product Datasheet

Radox - Integrin-Linked Kinase (ILK) (RCP9299)

Description	Integrin-linked kinase (ILK1) is a serine/threonine protein kinase, containing 4 ankyrin-like repeat...
Preservatives	25 mM Sodium acetate pH 4.8, 50 % Glycerol
Form/Appearance	Liquid
Storage	at -20 °C
Tag	N terminal containing a HIS tag
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	55.92 kDa
Uniprot ID	Q13418
Expiration Date	6 months from date of receipt.

