

## Product Datasheet

### Lipomaster 2000 Transfection Reagent (TL201-02)

|                          |                                                                                                                                                                                                                                                                                            |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Lipomaster 2000 Transfection Reagent is a liposome-mediated transfection reagent. It is applicable f...                                                                                                                                                                                    |
| <b>Storage</b>           | Store a 2 ~ 8°C. Adjust the shipping method according to the destination.                                                                                                                                                                                                                  |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                      |
| <b>Application notes</b> | Broad cell compatibility: Suitable for transfection of various adherent/suspension cells or co-transfection of siRNA and plasmid DNA<br>High transfection efficiency: Efficient transfection for both plasmid DNA and RNA<br>Low cytotoxicity: No need to change medium after transfection |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                                                            |

