

Product Datasheet

Lipomaster 3000 Transfection Reagent (TL301-02)

Description	Lipomaster 3000 Transfection Reagent is a liposome-mediated transfection reagent. It is applicable f...
Storage	Store a 2 ~ 8°C. Adjust the shipping method according to the destination.
Note	For research use only
Application notes	Superior performance: Efficient transfection of large plasmids and difficult-to-transfect cells Broad compatibility: Only one reagent can be used for transfection of plasmid DNA, mRNA, siRNA, and co-transfection, etc. Keep cell viability: Gentle on cells, low toxicity
Expiration Date	12 months from date of receipt.

