

Product Datasheet

Human NFKB1 ELISA Kit (dc654062-fb84-4232-9ac4-865f9b2d783a)

Description	Human NFKB1 ELISA Kit	Image not found or type unknown
Reactivity	Human	Human
Range	0.312-20ng/ml	NFKB1(Nuclear factor NF- kappa-B p1...
Target	NFKB1	
Storage	2-8 °C for 6 months	
Note	For research use only	
Uniprot ID	P19838	
Sensitivity	0.188ng/ml	
Expiration Date	Please enquire.	