



Product Datasheet

Nucleolar protein 3 antibody (orb100019)

Description	Nucleolar protein 3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse Nucleolar protein 3/Apoptosis repressor with CARD (1-100/220 aa)
Target	Nucleolar protein 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	O60936
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Flow cytometric analysis of A549 Cell suspension.