

Product Datasheet

2, 4-Dichlorophenoxyacetic acid antibody (orb10004)

Description	2, 4-Dichlorophenoxyacetic acid antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated 2,4-D
Target	2,4-D
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	94-75-7
Expiration Date	12 months from date of receipt.