



Product Datasheet

Insulin Receptor R antibody (orb100192)

Description	Insulin Receptor R antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Insulin receptor-related protein beta chain (801-900/1297 aa)
Target	Insulin Receptor R
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	129 kDa
Uniprot ID	P14616
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Flow cytometric analysis of A549 cells showing two populations, M1 and M2.