



Product Datasheet

HPGD antibody (orb100318)

Description	HPGD antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Prostaglandin dehydrogenase 1 (55-160/266 aa)
Target	Prostaglandin dehydrogenase 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	P15428
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse pl...



Immunohistochemical staining of Mouse pl...