



Product Datasheet

SQSTM1 antibody (orb100344)


Description	SQSTM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human SQSTM1/p62 (51-150/440 aa)
Target	SQSTM1/P62
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60 kDa
CAS Number	Sequestosome-1
Uniprot ID	Q13501
Dilution Range	IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of HepG2 cell using SQSTM1/...



IHC-P analysis of Mouse brain using SQST...



IHC-P analysis of rat stomach using SQST...