



Product Datasheet

AD7C-NTP antibody (orb10039)


Description	AD7C-NTP antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Neuronal thread protein AD7c-NTP (301-375/375 aa)
Target	Neuronal thread protein AD7c-NTP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1ug/Test, ELISA=1:5000-10000, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of mouse brain tis...



Immunohistochemical analysis of Mouse br...