



Product Datasheet

LIPG antibody (orb100394)


Description	LIPG antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LIPG. Please contact us for the exact immunogen sequence. The peptide is available as orb374809.
Target	LIPG. Recognized endogenous levels of total LIPG.
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	57 kDa
Uniprot ID	Q8VBX1, Q9Y5X9
Entrez	9388
Dilution Range	WB: 1:100-500, IHC-P: 1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human th...



Western blot analysis of rat kidney (lan...



IHC-P image of rat ovary tissue using an...