


Product Datasheet

ADAMTS1 antibody (orb10040)

Description	ADAMTS1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse ADAMTS1 (851-967/967 aa)
Target	ADAMTS1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	106 kDa
Uniprot ID	P97857
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Paraffin-embedded human lung carcinoma t...