



Product Datasheet

CRLF2 antibody (orb100404)

Description	CRLF2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse CRLF2 (251-350/370 aa)
Target	CRLF2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	Q8CII9
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



WB analysis of mouse BV2 cell (Lane1), m...



WB analysis of mouse kidney tissue...