

Product Datasheet

Clathrin heavy chain antibody (orb100509)

Description	Clathrin heavy chain antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Clathrin heavy chain. (401-500/1675 aa)
Target	Clathrin heavy chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	192 kDa
Uniprot ID	Q00610
Dilution Range	WB=1:500-2000, ICC=1:25
Expiration Date	12 months from date of receipt.