

Product Datasheet

Estrogen Receptor beta 2 antibody (orb100531)

Description

Estrogen Receptor beta 2 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Estrogen Receptor beta 2. Please contact us for the exact immunogen sequence. The peptide is available as orb12934.

Target

Estrogen Receptor beta 2

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

59 kDa

Uniprot ID

[P11475](#), [Q61539](#), [O95718](#)

NCBI

[004452](#), [004443](#)

Entrez



[2103](#)

Dilution Range


IF/ICC: 1:50-400, IHC-P: 1:50-400, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



IHC-P image of rat ovary tissue using Es...



IHC-P image of mouse testis tissue using...



Immunohistochemical staining of rat ovar...