


Product Datasheet


RNA polymerase II CTD repeat YSPTSPS (phospho-Ser5) antibody (orb100549)

Description

RNA polymerase II CTD repeat YSPTSPS



Immunohistochemical analysis of formalin...



Flow cytometric analysis of Blank contro...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC

Applications

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human RNA polymerase II CTD around the phosphorylation site of Ser5 (PT(p-S)PS)

Target

RNA polymerase II CTD repeat YSPTSPS (phospho S5)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

240 kDa

Uniprot ID

[P24928](#)

Dilution Range

Flow-Cyt=1ug/Test, WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.