



Product Datasheet

Syncytin 1 antibody (orb100573)


Description	Syncytin 1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Syncytin 1. Please contact us for the exact immunogen sequence. The peptide is available as orb374808.
Target	Syncytin 1
Form/Appearance	10mM PBS, 0.02 NaN ₃ , 2% rAlbumin
Concentration	100 ug (in 200 ul): 0.5 mg/ml 200 ug (in 400 ul): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	60 kDa
Uniprot ID	Q9UQF0
NCBI	001130925 , 001124397 , 11
Entrez	30816
Dilution Range	IHC-P: 1:200, WB: 1:200-800
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain (lane 1)



IHC-P staining of pig heart tissue using...



WB analysis of mouse brain (lane 1), rat...