



Product Datasheet

MSR1 antibody (orb100585)


Description	MSR1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	KLH conjugated synthetic peptide derived from human MSR1/CD204 (211-320/451 aa)
Target	MSR1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	P21757
Dilution Range	WB=1:500-2000, IF=1:100-500, Flow-Cyt=1ug/Test, IHC-P=1:100-500, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Paraffin-embedded human glioma tissue fi...



IHC-P of mouse lung tissue (MSR1 antibod...



Western blot analysis of human colon car...