

Product Datasheet

4 Hydroxynonenal antibody (orb100588)

Description	4 Hydroxynonenal antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated 4-Hydroxynonenal
Target	4 Hydroxynonenal
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	0.156 kDa
CAS Number	75899-68-2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.