

Product Datasheet

Pan 14-3-3 antibody (orb1006)

Description	Pan 14-3-3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human pan 14-3-3 (151-255/255 aa)
Target	14-3-3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	P63104
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1µg/Test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.