


Product Datasheet


TNFAIP5 antibody (orb100676)

Description	TNFAIP5 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human TNFAIP5/Pentraxin 3 (151-250/381 aa)
Target	TNFAIP5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P26022
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Paraffin-embedded rat lung tissue fixed

...



Western blot analysis of mouse embryo ly...