

Product Datasheet

Anti-Phospho-PAK6(Ser560)/APC (orb1006763)

Description

Anti-Phospho-PAK6(Ser560)/APC

Species/Host

Rabbit

Reactivity

Feline, Human, Mouse, Rat

Conjugation

APC

Tested

FC, ICC

Applications

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PAK6 around the phosphorylation site of Ser560

Target

PAK6

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

75kDa

Uniprot ID

Q9NQU5

Dilution Range

ICC=1:100

Expiration Date

12 months from date of receipt.