



Product Datasheet

F-Actin antibody (orb10068)

Description	F-Actin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic MAP peptide derived from human beta-Actin (1-50/374 aa)
Target	F-Actin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	P60709
Dilution Range	WB=1:1000-5000, IHC-P=1:100-500, Flow-Cyt=3µg/test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of rat skeletal mu...