



Product Datasheet

ESRRG antibody (orb100723)


Description	ESRRG antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human ESRRG/Estrogen Related Receptor gamma (61-160/458 aa)
Target	Estrogen Related Receptor gamma
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P62508
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



WB analysis of human endometrial cancer ...



Western blot analysis of human mammary f...



WB analysis of 293T cells (lane 1), 293T...